

RX400 Software - Release 1.6.3

Installation

Make sure that you have installed the 1.5.x release for support of the new update file format
If not please follow the instructions here: [Release 1.5.1](#)

Download and install the following file using the system updater to update to 1.6.3 from 1.5.x:

- [rx400_release_1.6.3.rx400update](#)

! The system automatically reboots after the installation of this update !

Versions included in this release:

- Drivecopy Service: v3.12.1
- Webinterface: v3.16.0
- Devicemanager: v1.7.0
- Updater Service: v2.0.0
- Worm Tools: v1.2.0
- Mediacreator: v1.1.2

Recovery Image

The current recovery image is version 1.5.1
You can find it here: [Release 1.5.1](#)

Direct download of current recovery image:

- [recovery_rx400medx_1.5.1.zip](#)

Changelog

RX400 Devicemanager (x5-devicemanager)

v1.7.0

- NEW: improved firewall

v1.6.0

- NEW: Support for Rimage Maestro Sticks on both variants (MedX & Central) [X5-604]

Drivecopy Service (x5-drivecopy)

v3.12.1

- FIX: Removed Rimage generic sticks (Product ID 6300) from Central support, now supported by generic license

v3.12.0

- FIX: Fix in issue that caused USB sticks to disappear and immediately reappear

v3.11.1

- FIX: Fix an issue which causes the touch screen to hang after long idle times
- FIX: Fix an issue which caused an error if a ticket that was still preparing was deleted

v3.11.0

- NEW: Support for Rimage Maestro Sticks on both variants (MedX & Central) [X5-604]
- NEW: Support for new Central usb media with new ids
- NEW: Support for sha512 hash algorithm for file verification
- FIX: Sometimes special filesystem files and directories caused verification errors
- FIX: In autoassign mode base line pro settings were set incorrect
- NEW: Extended logfile download

RX400 Webinterface (x5-web)

v3.16.0

- NEW: Support for RX400 DICOM extension

v3.15.0

- NEW: Extended simple ticket wizard [X5-607]
- NEW: Option to specify target directory name for read jobs [X5-603]
- FIX: There was an error with the simple wizard if you tried to create a ticket without label